LUMA



LUMA GO

User Manual

Scan to download the user manual

LUMA GO Portable Red Light Therapy Device

The LUMA GO is a portable and rechargeable red light therapy device, designed for convenient use at home or on the go.

Power Supply & Charging:

The device can be recharged using the provided USB cable. Simply connect it to any standard USB power source, such as the same USB outlet used for charging a mobile phone.

Battery Performance:

When fully charged, the LUMA GO provides approximately 2 hours of continuous operation, equivalent to 10–12 therapy sessions, depending on session length.


Dosage Recommendations

Distance: Hold the device approximately 4 to 12 inches (10–30 cm) away from the targeted area.
Frequency: It is recommended to perform red light therapy 2–3 times per week. For enhanced results, daily use is also suitable.
Multiple Sessions: If additional treatments are desired, the device may be used more than once per day. However, allow at least 6 hours between sessions to ensure your body fully processes the effects of the previous treatment.

Product Parameters




Wavelengths	660nm:850nm
Battery Size	7.4V 2500mAh
Battery Performance	2 hours
Dimensions	140 x 70 x 29 mm
Power Consumption	7.2W
Beam Angle	30 Degree
Timer	0-30 miniutes adjustable
Pulse	0-40 Hertz adjustable (only for NiR LEDs)
Charging time	1.5h (type-c)
Weight	244g
LED Lifespan	50, 000hrs



>200
58
40
30
Irradiance(mw/cm2)
1cm
5cm
10cm
15cm
Distance(cm)



After powering on, it displays for 1.5 seconds, then enters the main menu.



1. Press the **“Smart”** button to select Joint Care, Face & Eyes, Sleep, Skin, Workout, or Mood Health mode
2. Press the **“Custom”** button to choose RED, NIR, Timer, or Pulse.
3. Press **“ON”** to turn the mini panel on or off.




Press and hold **"Custom"** to turn mute on or off.



Select your preferred smart mode, then press "**ON**".

Press “ ▸ ” this icon to go to the next page.

Press this icon " ▬ " to return to the main screen.









1. After pressing **“Custom”** you can adjust the brightness settings for NIR or RED, as well as modify the time and pulse parameters.
2. Press “ ” this icon to select and access the settings for NIR, RED, time, and pulse.
3. Press this icon " " to return to the main screen.



Operating Screen



Pause Screen



End Screen



Since this is a handheld therapy light, a lock button prevents accidental touches. A quick press locks or unlocks the screen protection.




CE
FC
FDA
RoHS